



We have the solutions!



# pH indicator strips & papers



## pH strips

Product Code	Product Name	Graduation	Graduation in pH units	Pack size
CL09.1601	Universal Indicator Strips	pH 0 - 14	1	100 strips/pack
CL09.1602	pH Indicator Strips	pH 0 - 6.0	0.5	100 strips/pack
CL09.1603	pH Indicator Strips	pH 2.0 - 9.0	0.5	100 strips/pack
CL09.1604	pH Indicator Strips	pH 4.0 - 10.0	0.5	100 strips/pack
CL09.1605	pH Indicator Strips	pH 6.5 - 10.0	0.5	100 strips/pack
CL09.1606	pH Indicator Strips	pH 7.0 - 14.0	0.5	100 strips/pack
CL09.1607	pH Indicator Strips	pH 0 - 1.5	0.2/0/3	100 strips/pack
CL09.1608	pH Indicator Strips	pH 2.5 - 4.5	0.3	100 strips/pack
CL09.1609	pH Indicator Strips	pH 4.0 - 7.5	0.3	100 strips/pack

## pH papers

Product Code	Product Name	Graduation	Graduation in pH units	Pack size
CL09.1610	Universal Indicator paper	pH 1 - 14	1	5 m roll
CL09.1611	Universal Indicator paper	pH 1 - 14	1	3 × 5 m refill roll
CL09.1612	Universal Indicator paper	pH 1 - 11	1	5 m roll
CL09.1613	Universal Indicator paper	pH 1 - 11	1	3 × 5 m refill roll
CL09.1614	Comparator Paper	pH 0.5 - 5.5	0.5	5 m roll
CL09.1615	Comparator Paper	pH 0.5 - 5.5	0.5	3 × 5 m refill roll
CL09.1616	Comparator Paper	pH 5.5 - 9	0.5	5 m roll
CL09.1617	Comparator Paper	pH 5.5 - 9	0.5	3 × 5 m refill roll

FORR-LAB KFT.  
telephely  
1141 Budapest  
Fogarasi út 218

Tel. +36 1 273 1378  
Fax. +36 1 273 1379  
info@forrlab.hu  
www.forrlab.hu

FORR-LAB KFT.  
székhely  
1131 Budapest  
Madarász Viktor u. 1/A



We have the solutions!



# Chem-QUANT® Test strips



### Chem-QUANT® Test strips:

- Reliable and precise
- Non-bleeding
- Easy for on-site analysis
- Very quick

### Chem-QUANT® use:

- Dip the strip in the sample
- Wait 1 to 2 minutes for the reaction
- Compare with the color chart

Product code Chem-Lab	Product Name Chem-Lab	Graduation	Pack Size
CL09.0301	Chem-QUANT® Ammonium 400	0 - 10 - 25 - 50 - 100 - 200 - 400 mg/L NH4+	50 strips/pack
CL09.0302	Chem-QUANT® Carbonate Hardness	0 - 5 - 10 - 15 - 20 - 30 °d	100 strips/pack
CL09.0303	Chem-QUANT® Chlorine 20	0 - 1 - 3 - 5 - 20 mg/L	100 strips/pack
CL09.0304	Chem-QUANT® Chlorine 1000	0 - 50 - 100 - 250 - 500 - 1000 mg/L	100 strips/pack
CL09.0305	Chem-QUANT® Fluoride 100	0 - 10 - 25 - 50 - 100 mg/L	100 strips/pack
CL09.0306	Chem-QUANT® Glucose 2000	0 - 100 - 250 - 500 - 1000 - 2000 mg /dL	100 strips/pack
CL09.0307	Chem-QUANT® Nitrate 500	0 - 10 - 25 - 50 - 250 - 500 mg/L	100 strips/pack
CL09.0308	Chem-QUANT® Nitrite 80	0 - 2 - 5 - 10 - 20 - 40 - 80 mg/L	100 strips/pack
CL09.0309	Chem-QUANT® Peracetic Acid 50	0 - 5 - 10 - 20 - 30 - 50 mg/L	100 strips/pack
CL09.0310	Chem-QUANT® Peracetic Acid 2000	0 - 500 - 1000 - 1500 - 2000 mg/L	100 strips/pack
CL09.0311	Chem-QUANT® Peroxide 25	0 - 0.5 - 2 - 5 - 10 - 25 mg/L	100 strips/pack
CL09.0312	Chem-QUANT® Peroxide 1000	0 - 100 - 200 - 400 - 600 - 800 - 1000 mg/L	100 strips/pack
CL09.0313	Chem-QUANT® Phosphate 100	0 - 3 - 10 - 25 - 50 - 100 mg/L	50 strips/pack
CL09.0314	Chem-QUANT® QAC 25-1000	0 - 25 - 50 - 100 - 250 - 500 - 1000 mg/L Benzalkonium chloride	100 strips/pack
CL09.0315	Chem-QUANT® Silver 10	0 - 1 - 1.75 - 2.75 - 3.5 - 5 - 7 - 10 g/L	10 books of 20 strips
CL09.0316	Chem-QUANT® Sulphite 500	10 - 50 - 100 - 250 - 500 mg/L	100 strips/pack

FORR-LAB KFT.  
telephely  
1141 Budapest  
Fogarasi út 218

Tel. +36 1 273 1378  
Fax. +36 1 273 1379  
info@forrlab.hu  
www.forrlab.hu

FORR-LAB KFT.  
székhely  
1131 Budapest  
Madarász Viktor u. 1/A