

Púdermentes Nitril Gumikesztyű!

Elérhető a Forr-Lab Kft.-nél!!!



SATRA TECHNOLOGY Issued to: Master Glove Manufacturing Sdn Bhd
Lai Sui Tin
Jalan Haji Abdul Maman
8th Mile off Jalan Meru
41500 Klang
Selangor
Malaysia

Notified Body: 2777 SATRA customer number: P19025

EU Type-Examination Certificate

Certificate number: 2777/12716-01/E00-00

This EU Type-Examination Certificate covers the following product group(s) supported by testing to the relevant standards/technical specifications and examination of the technical file documentation:
Following the EU Type-Examination this product group has been shown to satisfy the applicable essential health and safety requirements of Annex II of the PPE Regulation (EU) 2016/425 as a Category III product.

Product reference: PFTN-FTB, PFTN-FTBP, PFTN-FTBK
Description: Five finger non sterile powder free nitrile examination glove, 3.2 mil

Classification:
EN ISO 374-1:2016 Type C Level EN 374-4: 2013
n-Heptane (J) 1 51.7%
Sodium hydroxide 40% (K) 5 -37.8%
Formaldehyde 37% (T) 5 -0.0%

EN ISO 374-5: 2016 Level
Protection against bacteria & fungi Pass
Protection against viruses Pass

Standards/Technical specifications applied:
EN ISO 374-1:2016, EN ISO 374-5:2016, EN 420: 2003+A1: 2009

Technical reports/Approval documents:
SATRA: CHM0277509/1844/JH/J, CHM0277509/1844/JH/J, CHM0277509/1844/JH/JAA/Final,
CHM0277509/1844/JPT/J
CENTEXBEL: 17.00609.02

Signed on behalf of SATRA: Tara Holly, Quincey Brown
Date first issued: 17/05/2019
Date of issue: 17/05/2019
Expiry date: 17/05/2024

Page 1 of 2
SATRA Technology Europe Limited, Brunner Business Park, Clonsilla, Dublin 15, Ireland

TECHNISCHES DATENBLATT **AMPri**

AMPri Handgloves, eBt 1 - Pfl. Urmacher-Allee 1a • 31425 Biele - Tel. 0541/7187-0 • Fax 0541/7187-17 • www.ampri.de • info@ampri.de

WHITE BASIC PLUS

Nitrile-Glove
powderfree / non-sterile
Article-number: 01200

Material: nitrile synthetical rubber
Powder: none
Surface: textured
Fit: ambidextrous
Colour: white
Elongation: before aging 6 N
after aging 6 N
Thickness: approx. 0,07 mm
Quality standard: ISO 9001:2000, ISO 13485:2003, CE
AQL 1.5, EN 420, EN 455, EN 374 1+2
Packing: 200 pcs per dispenser
10 dispenser per outer-carton

Sizes:	XS 5 - 5,5	S 6 - 6,5	M 7 - 7,5	L 8 - 8,5	XL 9 - 9,5
Width	75 mm ± 5	85 mm ± 5	95 mm ± 5	105 mm ± 5	115 mm ± 5
Length	240 mm ± 10	240 mm ± 10	240 mm ± 10	240 mm ± 10	240 mm ± 10

03/1.1

További információért keresse honlapunkat, www.forrlab.hu, vagy írjon nekünk!
info@forrlab.hu